


LED


Phototherapy

ULTRA HIGH
INTENSITY

L-15



L-9



TRUSTED
IRRADIANCE

BYRD TM **MEDITECH**
Innovative Technologies



LED Doublesurface

Irradiance $> 35\mu\text{W}/\text{cm}^2/\text{nm}$

- Wide Effective Area
- Long Life – > 20,000 hrs.
- No UV/IR Radiation
- No Heat Generation
- Low Power Consumption
- Negligible Maintenance Cost
- Light Weight & Compact
- Flexible Neck for Easy Use
- Height and Angle Adjustment Facility
- Mounted on Castor Wheels for Easy Mobility
- Can be used in Conjunction with any make Radiant Warmer
- Digital Time Totalizer
- Wave length : 400-550 nm
- Safety Cut Off in case of overheating (Optional)
- Total Hours and Therapeutic Hours Display
- 4 Castor Wheels
- White Light for Observation (Optional)